

NEW MULTIPLE ADJUSTMENT Patient Laser Goggles

Made of stainless steel construction with four adjustable features.
Includes white silicone eye cups and two retention straps.

Patent Pending

CE Approved



Nose Bridge Pivots Up and Down

Nose bridge pivots freely up and down without altering the position of the eye cups.



Eye Cups Tilt

Eye cups tilt up to 30 degrees providing a custom fit.




Removable Eye Cups with Silicone Eye Pads

Eye cups may be removed completely from the nose bridge to allow for individual treatment around the eyes.



Nose Bridge has Horizontal Adjustment

By pressing button, nose bridge easily adjusts horizontally to provide accurate pupillary distance setting for each patient.



KENTEK[®]

"Laser People Helping Laser People"[®]

Order part # KIP-CE

Call a Kentek Sales Representative: 1-800-432-2323

Secure Online Shopping at kenteklaserstore.com